


PIGTAIL
DRAINAGE
CATHETERS

BIOMETRIX
CRITICAL CARE SOLUTIONS



CLINICAL APPLICATIONS

Percutaneous Drainage of:

- Abscess and Fluid Collection
- Nephrostomy
- Biliary
- Pleural Effusion and Pneumothorax
- Pericardiocentesis




FEATURES AND BENEFITS

- **Radiopaque Thermoplastic Polyurethane**
Contains features that enable easy catheterization
and Traceable under reflection
- **Hydrophilic Coating**
Reduces coefficient friction
and facilitates easy catheter introduction

FEATURES AND BENEFITS


Pigtail Distal End

Prevents the catheter from sliding back



Safety String Lock Design

For strengthening fixation of pigtail




FEATURES AND BENEFITS

Large Oval Drainage Holes

Increase drainage efficiency


Holes Located Inside the Pigtail Loop

Prevents occlusion by tissue and reduces trauma




PLACEMENT TECHNIQUE

One Step Technique



PLACEMENT TECHNIQUE

Seldinger Technique



Guide wire

SETS AND ACCESSORIES

Biometrix offers customized trays and kits

